



MuDanJiang Plasma Physical Application Technology Co., Ltd.

No.99, ShengLiRing Road, JiangNan Development Area, Dong'An District,
MuDanJiang, HeiLongJiang, P.R. China Tel. +86-453-648-1685 Fax. +86-453-648-1766

Certificate of Analysis

Description: **CASP-120/120D**

Spec. \ Model	CASP-120	CASP-120D
Electronic Requirement	3P 380 \pm 38V 50/60HZ	
Power Consumption	4.0KW	
Overall Size	L845 x W1110 x H1780mm	L1100 x W1110 x H1780mm
Weight	470KG	515KG
Chamber Size	L682 x W611 x H362mm	L920 x W611 x H362mm
Chamber Volume	151 liters	203 liters
Operation Temperature	10 ~ 40 $^{\circ}$ C	
Operation Humidity	Under 85%(under less than 30 $^{\circ}$ C)	
Operation Pressure	70 ~ 106kPa	
Storage Temperature	-20 ~ 60 $^{\circ}$ C	
Storage Humidity	No special requirements, but as possible as keep dry	
Cycle Temperature Range	45 \pm 5 $^{\circ}$ C	
Sterilization Cycle Time	Fast: 22 min	Normal: 47 min Super: 57 min
Data Storage	60,000 cycles	

Goods are manufactured in compliance with our own manufacturing facility's quality standards and CE guidelines and all appropriate subsections.

Inspector : Han Wei/QC

Re-inspector : Zou Chen Ming/TEO

MuDanJiang Plasma Physical Application Technology Co., Ltd.



CASP CASP CASP CASP CASP CASP CASP CASP CASP CASP CASP