



MuDanJiang Plasma Physical Application Technology Co., Ltd.

No.99, ShengLiRing Road, JiangNan Development Area, Dong'An District,  
MuDanJiang, HeiLongJiang, P.R. China Tel. +86-453-648-1685 Fax. +86-453-648-1766

## Certificate of Analysis

Description: **CASP-50Flash**

Spec. \ Model	CASP-50Flash
Electronic Requirement	AC220 $\pm$ 22V 50/60HZ
Power Consumption	3.0KW
Overall Size	L681 x W612 x H794mm
Weight	138KG
Chamber Size	L551 x W227 x H210mm/U, L321 x W227 x H286mm/D
Chamber Volume	47 liters
Operation Temperature	10 ~ 40 $^{\circ}$ C
Operation Humidity	Under 85%(under less than 30 $^{\circ}$ C)
Operation Pressure	70 ~ 106kPa
Storage Temperature	-20 ~ 60 $^{\circ}$ C
Storage Humidity	No special requirements, but as possible as keep dry
Cycle Temperature Range	45 $\pm$ 5 $^{\circ}$ C
Sterilization Cycle Time	Fast: 18 min Normal: 42 min Super: 52 min
Data Storage	60,000 cycles

Goods are manufactured in compliance with our own manufacturing facility's quality standards and CE guidelines and all appropriate subsections.

Inspector : Han Wei/QC

Re-inspector : Zou Chen Ming/TEO

MuDanJiang Plasma Physical Application Technology Co., Ltd.



**CASP CASP CASP CASP CASP CASP CASP CASP CASP CASP CASP CASP**